



Fiber Optic Dual Switch with TTL Line

The FOS-2X2-TTL Fiber Optic Dual Switch was designed to provide you with flexibility in routing, splitting and controlling light. The FOS has two light channels. You can opt to have one light channel open at a time or have both closed. The FOS is useful for monitoring the drift of the light source or for measuring two samples with one spectrometer channel and one light source. The diagram below is an example of how the FOS can be utilized. In this setup, the FOS eliminates the need for a second spectrometer channel.

Item Code: FOS-2X2-TTL



In this setup, light enters a Bifurcated Optical Fiber Assembly and then splits into two arms, one for each sample. Light interacts with each sample and travels through another fiber assembly, each into its own port in the FOS. Another Bifurcated Assembly collects the light from the FOS and sends it to the spectrometer. Here you would switch the shutter on the FOS from one light channel to another in order to get clean data from each sample. Without the FOS, you would need another spectrometer channel to monitor the two samples.


Electronic TTL Shutter

In spectrometer setups, the INLINE-TTL TTL-driven shutter allows you to block the light path without disturbing the experiment -- for example, by turning the light source on and off. The laser-cut shutter is installed between two collimating lenses, which attach to two optical fibers. The INLINE-TTL is driven by a small board with a TTL input. Included is a cable for interfacing to a spectrometer.

Item Code: INLINE-TTL-S

Specifications

Dimensions:	140 mm x 50 mm x 50 mm
Weight:	~600 g
Shutter-Input:	TTL maximum 5 Hz
Power requirements:	12 VDC (power supply included)
Power consumption:	100 mA maximum
Maximum frequency:	5 Hz



Field of View Control

The Gershun Tube Kit (GER-KIT) controls the field of view of our SMA 905-terminated optical fiber. It also directly attaches to a spectrometer with an SMA 905 Connector. User-interchangeable apertures provide many different fields of view from 1° to 28°. (When the GER-KIT is used with our optical fiber, the field of view cannot exceed the optical fiber's 25° field of view if you are measuring radiance.)

Item Code: GER-KIT

Specifications

Material:	Black anodized aluminum
Interior:	Bead-blasted surface to reduce off-axis reflections
Connection:	Directly attaches to one of our spectrometers or couples to an SMA 905-terminated optical fiber with included adapter barrel
Apertures:	1°, 3°, 8°, 10° and 14° apertures included, providing 1°, 2°, 3°, 6°, 8°, 10°, 14°, 16°, 20° and 28° fields of view

