



# SpectroPipetter Microcell



## Easy to Use -- Just Pipette and Measure

The PIP-10-2 SpectroPipetter Microcell is a combination micropipetter and 10-mm pathlength microcell for low-volume sampling.


Samples are loaded into a capillary tube with an optical fiber plunger, which is activated by depressing the thumbpad and releasing it to draw in the fluid. A mirror on the distal side of the capillary completes the optical path.

## Requires 2 $\mu\text{L}$ of Sample

The pipetter is equipped with a bifurcated fiber, which couples to our spectrometers and compact light sources to create low-volume absorbance systems. The SpectroPipetter requires only 2  $\mu\text{L}$  of sample for a spectral measurement.

## Cleaning Kit


To remove fluid or dye from the cell, use the PIP-UCK Cleaning Kit. It contains an ultrasonic cleaning bath (at right) and a bottle of cleaning solution. Additional PIP-UCK-CS Cleaning Solution (below right) is also available.



### Specifications

|                                 |  |
|---------------------------------|--|
| Wavelength range:               | 230-850 nm                                     |
| Sample volume:                  | 2 $\mu\text{L}$                                |
| Light pathlength of cell:       | 10 mm  |
| Temperature range:              | 4 $^{\circ}\text{C}$ to 99 $^{\circ}\text{C}$  |
| Optical fiber length:           | 1.3 m  |
| Core diameter of launch fibers: | Bundle of (3) 200 $\mu\text{m}$ optical fibers |
| Core diameter of return fiber:  | 200 $\mu\text{m}$                              |

# Longpass Flow Cells



## 100x Increase in Sensitivity

LPC Longpass Flow Cells couple to our spectrometers and light sources for simple, efficient measurements of low-volume, low-concentration aqueous samples. With the LPC-1, you have a 1-meter cell with an internal volume of only 240  $\mu\text{L}$ , giving you 100x the sensitivity over a 1-cm pathlength cuvette holder as your sampling device.

## Easy to Use

LPC cells use a capillary tube as both the sample compartment and the light waveguide. You inject the sample into the fluidic ports with a syringe or pump; optical fibers connect to SMA 905 Connectors to deliver and return light to the spectrometer. We offer these cells in 1- and 5-meter pathlengths (call for other pathlengths). A 5-meter cell (250  $\mu\text{L}/\text{meter}$ ) increases the absorbance signal 500x more than a 1-cm cuvette. Also available is the LPC-CLEANKIT (see inset), a waveguide cleaning kit for the LPCs.

### Specifications

|                      |  |                             |   |
|----------------------|--|-----------------------------|---|
| Dimensions:          | 254 mm x 279 mm                            | Maximum sample temperature: | 160 $^{\circ}\text{C}$                          |
| Weight:              | 140 g                                      | Tubing inner diameter:      | 550 $\mu\text{m}$                               |
| Wavelength range:    | 230-800 nm for LPC-1, 325-700 nm for LPC-5 | Tubing:                     | Fused silica inner tubing coated with Teflon AF |
| Tubing volume:       | 250 $\mu\text{L}/\text{meter}$             | Fluid fittings:             | 1/16", 1/32" compression fittings               |
| Fiber connectors:    | SMA 905                                    | Maximum pressure:           | 2000 psi  |
| Fiber core diameter: | 400 $\mu\text{m}$                          | Chemical resistance:        | Most organic and inorganic solvents             |

